

Standard cleaning vs a disinfection cap for the decontamination of needle-free connectors.

AUTHORS:

Casey AL, Karpanen TJ, Nightingale P, Elliott TSJ.

University Hospitals Birmingham NHS Foundation Trust, Birmingham, UK.

Introduction

- A major factor attributed to the level of infection risk associated with needle-free connectors is the efficacy of disinfection of the injection ports [1].
- Disinfection of connectors may not intuitive which may lead to non-compliance [1].
- Caps which attach to injection ports of connectors incorporating disinfectants have been developed.
- These caps act as passive disinfection devices ensuring connectors are always clean.
- Clinical studies evaluating these devices demonstrated reductions in hub colonisation [2], and central-line associated bloodstream infections (CLABSI) [3-7].
- The aim of the study was to determine under controlled laboratory conditions whether a continuous passive disinfection cap containing 70% (v/v) IPA was as effective as defined cleaning with a 2% (w/v) CHG in 70% (v/v) IPA wipe.

Methods

- Connectors used: MicroClave™ neutral displacement connector (ICU Medical) and Care-Site™ positive-displacement connector (BBraun).
- Cleaning devices used: Curoso® caps containing 70% (v/v) IPA (3M Healthcare, figure) and 2% (w/v) CHG in 70% (v/v) IPA wipes (Sani-cloth CHG 2%, PDI).
- The injection port was inoculated with 5×10^6 CFU *Staphylococcus aureus* NCTC 6538 and allowed to dry.
- Disinfection caps were attached and compared with those cleaned with a wipe.
- All connectors were left at 20°C in air for 1, 3 or 7 days.
- Following contamination with *S. aureus*, a proportion of each type of connector were cleaned as above for 15 s with a wipe and allowed to dry for 30 s. These were then left for 7 days at 20°C and cleaned again with a wipe.
- Connectors were immersed into bijoux containing 1 mL of neutralizing solution
- The bijoux were sonicated for 10 min at 50 Hz and the solution inoculated onto chromogenic *S. aureus* plates (ChromID *S. aureus* [Biomerieux]).
- Median \log_{10} colony-forming unit (CFU) reductions and 95% confidence interval (CI) were calculated and data analyzed using the Mann-Whitney test. The level of significance was < 0.05 .

Curoso® cap containing 70% (v/v) IPA



Results

- The minimum CFU count on the control connectors (those which were not decontaminated after inoculation with *S. aureus*) was $5.17 \log_{10}$ CFU for MicroClave™ and $5.49 \log_{10}$ CFU for CareSite™.
- Total kill (TK) therefore always represented a ≥ 5.17 or $\geq 5.49 \log_{10}$ CFU reduction, respectively.
- \log_{10} CFU reductions in *S. aureus* following decontamination for 15 s with a wipe and incubation at room temperature for 1, 3 or 7 days or application of the disinfection cap for 1, 3 or 7 days is shown in the table.
- The disinfection cap resulted in a significantly higher reductions in *S. aureus* than the wipe.
- There was no difference in the \log_{10} CFU reduction of *S. aureus* between the two different types of connector.
- Decontamination of the MicroClave™ with a wipe following inoculation with *S. aureus* and following each subsequent incubation period resulted in a higher \log_{10} CFU reduction as compared to cleaning only following contamination ($p = 0.009$).
- The disinfection cap resulted in a significantly higher \log_{10} CFU reduction compared to the two decontaminations with wipes for both MicroClave™ ($p = 0.041$), and CareSite™ ($p < 0.0001$).

Median (95%CI) \log_{10} *Staphylococcus aureus* CFU reductions on two types of needle-free connectors after 1, 3 and 7 days following two decontamination methods. Total kill (TK) = no growth of *S. aureus* following decontamination.

Day	Decontamination method	Connector studied: MicroClave®	Comparison of wipe vs disinfection cap (P value)	Connector studied: CareSite®	Comparison of wipe vs disinfection cap (P value)	Comparison of MicroClave® vs CareSite® (P value)
1.	2% (w/v) CHG in 70% (v/v) IPA wipe	>6.45* (4.97-TK)	<0.0001*	TK (4.29-TK)	<0.0001*	0.49
	Disinfection cap	TK (TK-TK)		TK (TK-TK)		0.75
3.	2% (w/v) CHG in 70% (v/v) IPA wipe	4.66 (4.34-4.95)	<0.0001*	4.77 (4.39- 5.68)	<0.0001*	0.98
	Disinfection cap	TK (TK-TK)		TK (TK-TK)		0.057
7.	2% (w/v) CHG in 70% (v/v) IPA wipe	TK (TK-TK)	<0.0001*	TK (5.20-TK)	<0.0001*	0.15
	Disinfection cap	TK (TK-TK)		TK (TK-TK)		1.00

*=the median was half-way between the values of 6.45 and TK. +The reductions were greater for the disinfection cap.

Discussion & Conclusion

- Under controlled laboratory conditions a disinfection cap containing 70% (v/v) IPA was more effective at reducing microbial contamination of injection ports of needle-free connectors when compared to cleaning with 2% (w/v) CHG in 70% (v/v) IPA wipes.
- This may reflect the continuous antimicrobial activity of the decontamination offered by the caps.
- It could be argued that given the significant \log_{10} CFU reductions observed with the wipe in this study, there is no requirement for the disinfection cap. However, if wipe compliance is low, the disinfection caps could prove a useful tool.
- The results of this study support the SHEA/IDSA practice special approach recommendation for preventing CLABSI to 'use an antiseptic-containing hub/connector cap/port protector to cover connectors' [8].

Acknowledgements

- We would like to thank Karen Burgess for her assistance in the laboratory.
- Funding for this study was provided by 3 M.

References

1. Moureau NL, Flynn J. Disinfection of needleless connector hubs: clinical evidence systematic review. *Nurs Res Pract* 2015;article ID 796762.
2. Gutiérrez Nicolás F, Nazco Casanigo GJ, Viña Romero MM, Gonzalaz Garcia J, Ramos Diaz R, Perez Perez JA. Reducing the degree of colonisation of venous access catheters by continuous passive disinfection. *Eur J Hosp Pharm*. 2016;23:131-3
3. Merrill KC, Sumner S, Linford L, Taylor C, Macintosh C. Impact of universal disinfectant cap implementation on central-line associated bloodstream infections. *Am J Infect Control*. 2014;40(12):1274-7.
4. Sweet MA, Cumpston A, Briggs F, Craig M, Hamadani M. Impact of alcohol impregnated port protectors and needleless neutral pressure connectors on central line associated bloodstream infections and contamination of blood cultures in an inpatient oncology unit. *Am J Infect Control*. 2012;40(10):931-4.
5. Ramirez C, Lee AM, Welch K. Central venous catheter protective connector caps reduce intraluminal catheter-related infection. *J Assoc Vasc Access*. 2012;17(4):210-3.
6. Stango C, Ruryan D, Stern J, Macri I, Vacca M. A successful approach to reducing bloodstream infections based on a disinfection device for intravenous needleless connector hubs. *J Infus Nurs*. 2014;37(6):462-5.
7. Wright MO, Tropp J, Schora DM, et al. Continuous passive disinfection of catheter hubs prevents contamination and bloodstream infection. *Am J Infect Control*. 2013;41:33-8.
8. Marshall J, Mermel LA, Fakih M, et al. Strategies to prevent central line-associated bloodstream infections in acute care hospitals: 2014 update. *Infect Control Hosp Epidemiol*. 2014;35(7):753-71.